


# 4°C Blood Bank Refrigerator

*Super High Quality*

*Precise Temperature Control*

*Safeguarding Security of Blood*



## XC - 268L



### Constant Temperature under Intelligent Control

- Return air duct design, ensuring a temperature control accuracy at 0.1°C inside the cabinet;
- High-precision computerized temperature control system, built-in sensors for upper/lower zone temperature, ambient temperature, evaporator temperature, and operation control, ensuring safe and stable operation.



### Security System

Perfect audible and visual alarm system comes with alarm functions for high temperature, low temperature, sensor fault, door opening, and power failure.



### High-efficiency Refrigeration

- Advanced air cooling design, precision temperature control, safeguarding security of blood;
- Stainless steel inner container, copper tube evaporator, powerful refrigeration.






### Precise Temperature Control

1-inch high-brightness digital temperature display, allowing to set the temperature precision within 2~6°C inside the cabinet, and the precision of temperature display accuracy at 0.1°C.



### Refrigeration System

- Equipped with brand high-efficient compressor, EBM fan motor, performing efficient energy and quiet.
- High-precision computerized temperature control system, built-in sensors for upper/lower zone temperature, ambient temperature, evaporator temperature, and operation control, ensuring safe and stable operation.



### DC/DC Power Management

- Equipped with DC/DC power management module, optimized energy saving;
- High-precision voltage control, stabilizing the voltage and refrigerator performance;
- Reduce the redundancy heat, completely eradicate overheating risk.



### Automatic Defrost

- Return air duct design, ensuring a temperature control precision within 2~6°C inside the cabinet;
- Equipped with automatic defrost function, force defrost.



### Alarm System

- Perfect audible and visual alarm system comes with alarm function for high temperature, low temperature, sensor failure, door ajar, and power failure.




### Human-oriented

- Equipped with 5 dip-coating shelves with label strips;
- 15 dip-coating blood bag frames (optional stainless steel frame), can hold 165 blood bags in 450ml for each.
- There are two acrylic inner doors, which help to better insulate and keep out bacteria.


## ✓ Scope of Application

Suitable for use in blood banks, hospitals, health and disease prevention systems, research institutions, colleges & universities, the electronic industry, biological engineering, laboratories in colleges & universities, military enterprises, deep-sea fishing companies, etc.

## 🔗 Performance Data / Cooling Curve



## 📏 External Dimensions



## Specification Chart

4°C Blood Bank Refrigerator	
Model	XC-268L
Cabinet Type	Upright
Capacity(L)	268
Internal Size(W*D*H)mm	530*490*1145
External Size(W*D*H)mm	640*700*1856
Package Size(W*D*H)mm	700*760*2013
NW/GW(Kgs)	125/161
Performance	
Temperature Setting	2~6°C
Ambient Temperature	16~32°C
Climate Class	N
Controller	Microprocessor
Display	Digital display
Refrigeration	
Compressor	1pc
Cooling Method	Air cooling
Defrost Mode	Automatic
Refrigerant	R134a
Insulation Thickness(mm)	R/L/U/B:54,D:55
Construction	
External Material	Cold rolled steel sheet
Inner Material	Stainless steel
Shelves	5(coated steel wired shelf)
Door Lock with key	Yes
Blood Basket	15
Access Port	1pc. Ø 25 mm
Casters	4 (2 casters with brake) + (2 leveling feet)
Inner door	2
Data Logging/Interval/Recording Time	USB/Record every 10 minutes / 2 years
Backup Battery	Yes
Alarm	
Temperature	High/Low temperature,High ambient temperature
Electrical	Power failure , Low battery
System	Sensor failure, Door ajar, Condenser overheating,Built-in USB datalogger failure,Communication error
Electrical	
Power Supply(V/HZ)	220-230/50
Bai VØ1 fi	400
Bai er Consumpf[a` /=-I Z!\$&Zfi	5.3
DSfed CgdW f/3fi	2.5
Accessories	
Standard	RS485, Remote alarm contact

\*The model, parameters and performance specified in this brochure may be changed without prior notice because of product upgrading.

\*There may be differences between the product images shown in this brochure and the actual products. When you are buying any product, please check the actual product.